



MAGEC3 Antibody, HRP conjugated

Product Code	CSB-PA819466LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD91
Immunogen	Recombinant Human Melanoma-associated antigen C3 protein (75-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen C3 (Cancer/testis antigen 7.2) (CT7.2) (Hepatocellular carcinoma-associated antigen 2) (MAGE-C3 antigen), MAGEC3, HCA2
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MAGEC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.