



MAOA Antibody, Biotin conjugated

Product Code	CSB-PA013398LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21397
Immunogen	Recombinant Human Amine oxidase [flavin-containing] A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amine oxidase [flavin-containing] A (EC 1.4.3.4) (Monoamine oxidase type A) (MAO-A), MAOA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MAOA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.