



MAPK13 Antibody, Biotin conjugated

Product Code	CSB-PA013452LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15264
Immunogen	Recombinant Human Mitogen-activated protein kinase 13 protein (1-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 13 (MAP kinase 13) (MAPK 13) (EC 2.7.11.24) (Mitogen-activated protein kinase p38 delta) (MAP kinase p38 delta) (Stress-activated protein kinase 4), MAPK13, PRKM13 SAPK4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAPK13