



MAPK13 Antibody, HRP conjugated

Product Code	CSB-PA013452LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15264
Immunogen	Recombinant Human Mitogen-activated protein kinase 13 protein (1-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 13 (MAP kinase 13) (MAPK 13) (EC 2.7.11.24) (Mitogen-activated protein kinase p38 delta) (MAP kinase p38 delta) (Stress-activated protein kinase 4), MAPK13, PRKM13 SAPK4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAPK13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.