



MAS1L Antibody, HRP conjugated

Product Code	CSB-PA013506LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35410
Immunogen	Recombinant Human Mas-related G-protein coupled receptor MRG protein (318-378AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mas-related G-protein coupled receptor MRG (MAS-R) (MAS1-like), MAS1L, MRG
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MAS1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.