




MED7 Antibody

Product Code	CSB-PA013667LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43513
Immunogen	Recombinant Human Mediator of RNA polymerase II transcription subunit 7 protein (1-233AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 7 (hMED7) (Activator-recruited cofactor 34 kDa component) (ARC34) (Cofactor required for Sp1 transcriptional activation subunit 9) (CRSP complex subunit 9) (Mediator complex subunit 7) (RNA polymerase transcriptional regulation mediator subunit 7 homolog) (Transcriptional coactivator CRSP33), MED7, ARC34 CRSP9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MED7

Image



Western blot

All lanes: MED7 antibody at 6µg/ml

Lane 1: Jurkat whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: Mouse thymus tissue

Lane 4: NIH/3T3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 28 kDa

Observed band size: 28 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.