



# MED7 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA013667LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43513
<b>Immunogen</b>	Recombinant Human Mediator of RNA polymerase II transcription subunit 7 protein (1-233AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Mediator of RNA polymerase II transcription subunit 7 (hMED7) (Activator-recruited cofactor 34 kDa component) (ARC34) (Cofactor required for Sp1 transcriptional activation subunit 9) (CRSP complex subunit 9) (Mediator complex subunit 7) (RNA polymerase transcriptional regulation mediator subunit 7 homolog) (Transcriptional coactivator CRSP33), MED7, ARC34 CRSP9
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	MED7
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.