



MFAP3L Antibody, Biotin conjugated

Product Code	CSB-PA013748LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75121
Immunogen	Recombinant Human Microfibrillar-associated protein 3-like protein (29-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microfibrillar-associated protein 3-like (Testis development protein NYD-SP9), MFAP3L, KIAA0626
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MFAP3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.