



MFSD5 Antibody, Biotin conjugated

Product Code	CSB-PA747224LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6N075
Immunogen	Recombinant Human Molybdate-anion transporter protein (216-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Molybdate-anion transporter (Major facilitator superfamily domain-containing protein 5) (Molybdate transporter 2 homolog) (hsMOT2), MFSD5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MFSD5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.