





# MICAL2 Antibody

<b>Product Code</b>	CSB-PA013808LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O94851
<b>Immunogen</b>	Recombinant Human [F-actin]-monooxygenase MICAL2 protein (1-495AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	[F-actin]-monooxygenase MICAL2 (EC 1.14.13.225) (Molecule interacting with CasL protein 2) (MICAL-2), MICAL2, KIAA0750 MICAL2PV1 MICAL2PV2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	MICAL2

## Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA013808LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA013808LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.