



# MORC2 Antibody, Biotin conjugated

Product Code	CSB-PA896555LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6X9
Immunogen	Recombinant Human MORC family CW-type zinc finger protein 2 protein (2-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORC family CW-type zinc finger protein 2 (Zinc finger CW-type coiled-coil domain protein 1), MORC2, KIAA0852 ZCWCC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MORC2