



MPPED1 Antibody, Biotin conjugated

Product Code	CSB-PA014767LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15442
Immunogen	Recombinant Human Metallophosphoesterase domain-containing protein 1 protein (1-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallophosphoesterase domain-containing protein 1 (EC 3.1.-.-) (Adult brain protein 239) (239AB), MPPED1, C22orf1 FAM1A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MPPED1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.