



MRC1 Antibody, Biotin conjugated

Product Code	CSB-PA014783LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22897
Immunogen	Recombinant Human Macrophage mannose receptor 1 protein (655-1213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage mannose receptor 1 (MMR) (C-type lectin domain family 13 member D) (C-type lectin domain family 13 member D-like) (Human mannose receptor) (hMR) (Macrophage mannose receptor 1-like protein 1) (CD antigen CD206), MRC1, CLEC13D CLEC13DL MRC1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MRC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.