



MRPL18 Antibody, FITC conjugated

Product Code	CSB-PA861987LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0U6
Immunogen	Recombinant Human 39S ribosomal protein L18, mitochondrial protein (14-180AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L18, mitochondrial (L18mt) (MRP-L18) (Mitochondrial large ribosomal subunit protein uL18m), MRPL18
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPL18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.