



MRPL22 Antibody, FITC conjugated

Product Code	CSB-PA878907LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWU5
Immunogen	Recombinant Human 39S ribosomal protein L22, mitochondrial protein (41-206AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L22, mitochondrial (L22mt) (MRP-L22) (39S ribosomal protein L25, mitochondrial) (L25mt) (MRP-L25) (Mitochondrial large ribosomal subunit protein uL22m), MRPL22, MRPL25 RPML25
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPL22
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.