



MRPL3 Antibody, Biotin conjugated

Product Code	CSB-PA014836LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09001
Immunogen	Recombinant Human 39S ribosomal protein L3, mitochondrial protein (201-348AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L3, mitochondrial (L3mt) (MRP-L3) (Mitochondrial large ribosomal subunit protein uL3m), MRPL3, MRL3 RPML3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPL3