




MRPL38 Antibody

Product Code	CSB-PA856920LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DV4
Immunogen	Recombinant Human 39S ribosomal protein L38, mitochondrial protein (35-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L38, mitochondrial (L38mt) (MRP-L38) (Mitochondrial large ribosomal subunit protein mL38), MRPL38
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPL38

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA856920LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.