



MRPL48 Antibody, FITC conjugated

Product Code	CSB-PA836217LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GC5
Immunogen	Recombinant Human 39S ribosomal protein L48, mitochondrial protein (29-212AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L48, mitochondrial (L48mt) (MRP-L48) (Mitochondrial large ribosomal subunit protein mL48), MRPL48
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPL48
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.