



MRPL53 Antibody, FITC conjugated

Product Code	CSB-PA853410LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EL3
Immunogen	Recombinant Human 39S ribosomal protein L53, mitochondrial protein (8-113AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L53, mitochondrial (L53mt) (MRP-L53) (Mitochondrial large ribosomal subunit protein mL53), MRPL53
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MRPL53
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.