



MRPS10 Antibody, Biotin conjugated

Product Code	CSB-PA014875LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82664
Immunogen	Recombinant Human 28S ribosomal protein S10, mitochondrial protein (1-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S10, mitochondrial (MRP-S10) (S10mt) (Mitochondrial small ribosomal subunit protein uS10m), MRPS10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MRPS10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.