



MRPS14 Antibody, FITC conjugated

Product Code	CSB-PA014882LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60783
Immunogen	Recombinant Human 28S ribosomal protein S14, mitochondrial protein (1-128AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S14, mitochondrial (MRP-S14) (S14mt) (Mitochondrial small ribosomal subunit protein uS14m), MRPS14
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPS14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.