




MRPS2 Antibody

Product Code	CSB-PA897102LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y399
Immunogen	Recombinant Human 28S ribosomal protein S2, mitochondrial protein (1-296AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S2, mitochondrial (MRP-S2) (S2mt) (Mitochondrial small ribosomal subunit protein uS2m), MRPS2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPS2

Image



Western blot

All lanes: MRPS2 antibody at 2.5µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: U251 whole cell lysate

Lane 3: MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.