



MRPS28 Antibody, HRP conjugated

Product Code	CSB-PA896495LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2Q9
Immunogen	Recombinant Human 28S ribosomal protein S28, mitochondrial protein (72-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S28, mitochondrial (MRP-S28) (S28mt) (28S ribosomal protein S35, mitochondrial) (MRP-S35) (S35mt) (Mitochondrial small ribosomal subunit protein bS1m), MRPS28, MRPS35
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPS28
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.