



MRPS31 Antibody, HRP conjugated

Product Code	CSB-PA846089LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92665
Immunogen	Recombinant Human 28S ribosomal protein S31, mitochondrial protein (66-395AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S31, mitochondrial (MRP-S31) (S31mt) (Imogen 38) (Mitochondrial small ribosomal subunit protein mS31), MRPS31, IMOGN38
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPS31
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.