



MRPS36 Antibody, Biotin conjugated

Product Code	CSB-PA014927LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82909
Immunogen	Recombinant Human 28S ribosomal protein S36, mitochondrial protein (1-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	28S ribosomal protein S36, mitochondrial (MRP-S36) (S36mt), MRPS36
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPS36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.