




MRRF Antibody

Product Code	CSB-PA850282LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96E11
Immunogen	Recombinant Human Ribosome-recycling factor, mitochondrial protein (56-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-recycling factor, mitochondrial (RRF) (Ribosome-releasing factor, mitochondrial), MRRF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRRF

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA850282LA01HU at dilution of 1:100