




MS4A3 Antibody

Product Code	CSB-PA856947LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HJ5
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 3 protein (1-49AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 3 (CD20 antigen-like protein) (Hematopoietic-specific transmembrane protein 4) (HTm4), MS4A3, CD20L HTM4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MS4A3

Image


Immunofluorescent analysis of U251 cells using CSB-PA856947LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.