



MS4A7 Antibody, Biotin conjugated

Product Code	CSB-PA875636LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZW8
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 7 protein (1-47AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 7 (CD20 antigen-like 4) (CD20/FC-epsilon-RI-beta family member 4) (Four-span transmembrane protein 2), MS4A7, 4SPAN2 CD20L4 CFFM4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MS4A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.