



MT1H Antibody, Biotin conjugated

Product Code	CSB-PA17399D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80294
Immunogen	Recombinant Human Metallothionein-1H protein (1-59AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallothionein-1H (MT-1H) (Metallothionein-0) (MT-0) (Metallothionein-IH) (MT-IH), MT1H
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MT1H
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.