



MT1M Antibody, Biotin conjugated

Product Code	CSB-PA17429D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N339
Immunogen	Recombinant Human Metallothionein-1M protein (1-61AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallothionein-1M (MT-1M) (Metallothionein-IM) (MT-IM), MT1M, MT1K
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MT1M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.