



MTHFR Antibody

Product Code	CSB-PA015158YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42898
Immunogen	Recombinant Human Methylenetetrahydrofolate reductase protein (1-656AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylenetetrahydrofolate reductase (EC 1.5.1.20), MTHFR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTHFR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.