





MXD3 Antibody

Product Code	CSB-PA015252LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW11
Immunogen	Recombinant Human Max dimerization protein 3 protein (1-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Max dimerization protein 3 (Max dimerizer 3) (Class C basic helix-loop-helix protein 13) (bHLHc13) (Max-associated protein 3) (Max-interacting transcriptional repressor MAD3) (Myx), MXD3, BHLHC13 MAD3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MXD3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA015252LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA015252LA01HU at dilution of 1:100