





# MXD3 Antibody

<b>Product Code</b>	CSB-PA015252LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BW11
<b>Immunogen</b>	Recombinant Human Max dimerization protein 3 protein (1-206AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Max dimerization protein 3 (Max dimerizer 3) (Class C basic helix-loop-helix protein 13) (bHLHc13) (Max-associated protein 3) (Max-interacting transcriptional repressor MAD3) (Myx), MXD3, BHLHC13 MAD3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	MXD3

## Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA015252LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA015252LA01HU at dilution of 1:100



---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.