




MYH7B Antibody

Product Code	CSB-PA015301LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7E2Y1
Immunogen	Recombinant Human Myosin-7B protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-7B (Antigen MLAA-21) (Myosin cardiac muscle beta chain) (Myosin heavy chain 7B, cardiac muscle beta isoform) (Slow A MYH14), MYH7B, KIAA1512
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH7B

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA015301LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA015301LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.