



MYH7B Antibody, Biotin conjugated

Product Code	CSB-PA015301LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7E2Y1
Immunogen	Recombinant Human Myosin-7B protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-7B (Antigen MLAA-21) (Myosin cardiac muscle beta chain) (Myosin heavy chain 7B, cardiac muscle beta isoform) (Slow A MYH14), MYH7B, KIAA1512
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH7B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.