



NANOG Antibody, HRP conjugated

Product Code	CSB-PA888008LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9S0
Immunogen	Recombinant Human Homeobox protein NANOG protein (1-305AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein NANOG (Homeobox transcription factor Nanog) (hNanog), NANOG
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	NANOG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.