



NAPA Antibody, HRP conjugated

Product Code	CSB-PA015447LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54920
Immunogen	Recombinant Human Alpha-soluble NSF attachment protein (1-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-soluble NSF attachment protein (SNAP-alpha) (N-ethylmaleimide-sensitive factor attachment protein alpha), NAPA, SNAPA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAPA