



NARS2 Antibody, Biotin conjugated

Product Code	CSB-PA856950LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96I59
Immunogen	Recombinant Human Probable asparagine--tRNA ligase, mitochondrial protein (15-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable asparagine--tRNA ligase, mitochondrial (EC 6.1.1.22) (Asparaginyl-tRNA synthetase) (AsnRS), NARS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NARS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.