



# NARS2 Antibody, HRP conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA856950LB01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q96I59   |
| <b>Immunogen</b>           | Recombinant Human Probable asparagine--tRNA ligase, mitochondrial protein (15-210AA)                       |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | HRP  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                           |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Probable asparagine--tRNA ligase, mitochondrial (EC 6.1.1.22) (Asparaginyl-tRNA synthetase) (AsnRS), NARS2 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Signal Transduction  |
| <b>Target Names</b>        | NARS2  |