



NAT8 Antibody, HRP conjugated

Product Code	CSB-PA887033LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHE5
Immunogen	Recombinant Human N-acetyltransferase 8 protein (64-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetyltransferase 8 (EC 2.3.1.-) (Acetyltransferase 2) (ATase2) (Camello-like protein 1) (CysteinyI-conjugate N-acetyltransferase) (CCNAT) (EC 2.3.1.80), NAT8, CML1 GLA TSC501
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAT8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.