



NDFIP1 Antibody, Biotin conjugated

Product Code	CSB-PA887119LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT67
Immunogen	Recombinant Human NEDD4 family-interacting protein 1 protein (2-116AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NEDD4 family-interacting protein 1 (Breast cancer-associated protein SGA-1M) (NEDD4 WW domain-binding protein 5) (Putative MAPK-activating protein PM13) (Putative NF-kappa-B-activating protein 164) (Putative NFKB and MAPK-activating protein), NDFIP1, N4WBP5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NDFIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.