



NIF3L1 Antibody, Biotin conjugated

Product Code	CSB-PA880926LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZT8
Immunogen	Recombinant Human NIF3-like protein 1 protein (1-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NIF3-like protein 1 (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 1 protein), NIF3L1, ALS2CR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NIF3L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.