



NKG7 Antibody, HRP conjugated

Product Code	CSB-PA613690LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16617
Immunogen	Recombinant Human Protein NKG7 protein (30-60AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NKG7 (G-CSF-induced gene 1 protein) (GIG-1 protein) (Granule membrane protein of 17 kDa) (GMP-17) (Natural killer cell protein 7) (p15-TIA-1), NKG7, GIG1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NKG7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.