



NKIRAS1 Antibody, Biotin conjugated

Product Code	CSB-PA885767LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYS0
Immunogen	Recombinant Human NF-kappa-B inhibitor-interacting Ras-like protein 1 protein (1-192AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NF-kappa-B inhibitor-interacting Ras-like protein 1 (I-kappa-B-interacting Ras-like protein 1) (Kappa B-Ras protein 1) (KappaB-Ras1), NKIRAS1, KBRAS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NKIRAS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.