



NOSTRIN Antibody, Biotin conjugated

Product Code	CSB-PA818236LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVI9
Immunogen	Recombinant Human Nostrin protein (1-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nostrin (BM247 homolog) (Nitric oxide synthase traffic inducer) (Nitric oxide synthase trafficker) (eNOS-trafficking inducer), NOSTRIN
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOSTRIN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.