





NPC2 Antibody

Product Code	CSB-PA015976LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61916
Immunogen	Recombinant Human NPC intracellular cholesterol transporter 2 protein (20-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NPC intracellular cholesterol transporter 2 (Epididymal secretory protein E1) (Human epididymis-specific protein 1) (He1) (Niemann-Pick disease type C2 protein), NPC2, HE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NPC2

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA015976LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA015976LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.