



NPC2 Antibody, Biotin conjugated

Product Code	CSB-PA015976LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61916
Immunogen	Recombinant Human NPC intracellular cholesterol transporter 2 protein (20-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NPC intracellular cholesterol transporter 2 (Epididymal secretory protein E1) (Human epididymis-specific protein 1) (He1) (Niemann-Pick disease type C2 protein), NPC2, HE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NPC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.