



NPC2 Antibody, HRP conjugated

Product Code	CSB-PA015976LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61916
Immunogen	Recombinant Human NPC intracellular cholesterol transporter 2 protein (20-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NPC intracellular cholesterol transporter 2 (Epididymal secretory protein E1) (Human epididymis-specific protein 1) (He1) (Niemann-Pick disease type C2 protein), NPC2, HE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NPC2