



NR0B2 Antibody, HRP conjugated

Product Code	CSB-PA624025LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15466
Immunogen	Recombinant Human Nuclear receptor subfamily 0 group B member 2 protein (1-257AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 0 group B member 2 (Orphan nuclear receptor SHP) (Small heterodimer partner), NR0B2, SHP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NR0B2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.