



NR0B2 Antibody, HRP conjugated

Product Code	CSB-PA624025LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15466
Immunogen	Recombinant Human Nuclear receptor subfamily 0 group B member 2 protein (1-257AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 0 group B member 2 (Orphan nuclear receptor SHP) (Small heterodimer partner), NR0B2, SHP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NR0B2